



Cyrusbioscience

Specification

Ponceau S Solution

0.2% (W/V) in 3% acetic acid

Product Code:	101-6226-79-5
CAS number:	6226-79-5
Formula:	C ₂₂ H ₁₂ N ₄ Na ₄ O ₁₃ S ₄
Molecular mass.:	760.6
Storage:	Store at Room Temperature

Specification

Appearance(color):	Dark red
Appearance(form):	Liquid
Suitability (Suitable for use in cellulose acetate electrophoresis):	Suitable
Solubility:	Soluble in water (1mg/ml). Insoluble in ether.

Responsible Laboratory Manager Quality Control: Taylor Lin

This document has been produced electronically and is valid without a signature.