



**Cyrusbioscience**

---

## Water Bath Cleaner®

**Cat. Number:** 2010

**Size:** 250mL

**Storage:**

Store at room temperature.

**Description**

Water Bath Cleaner® can keeps water baths, incubators, temperature controlled instruments, humidifiers, etc. clean and clear.

**Application**

Used to keep any circulating, temperature-controlled instrument clean, instruments like circulating Water baths, incubators, temperature controlled instruments, humidifiers, automatic analyzer baths and Blood/gas analyzers.

**Usage**

250mL of Water Bath Cleaner® can be used to treat about 1600L of water.

For 1L of water, add 0.15-0.2mL of Water Bath Cleaner®.

**Note! Not recommended for use with soap or anionic wetting agents.**