

堅持給您最好的 緩衝液原料&溶液



塞魯士生技緩衝液三大保證:

1. 高品質的原料
2. 嚴格的檢驗系統
3. 客製化服務
(特殊緩衝溶液配製及諮詢)

CAS NO.	品項	包裝	價格
7365-82-4	ACES	25g/ 100g	1,580/ 5,250
6976-37-0	Bis-Tris	100g/ 500g	2,100/ 8,400
64431-96-5	Bis-Tris propane	25g/ 100g	2,150/ 6,500
10191-18-1	BES	100g/ 500g	1,890/ 7,560
150-25-4	Bicine	100g/ 500g	1,050/ 4,200
1135-40-6	CAPS	100g/ 500g	1,350/ 5,500
75621-03-3	CHAPS	1g/ 5g	840/ 2,730
103-47-9	CHES	25g/ 100g	1,200/ 3,950
7365-45-9	HEPES	100g/ 500g	882/ 3,150
16052-06-5	HEPPS	25g/ 100g	1,470/ 4,200

好物推薦

HEPES緩衝液—

細胞培養中最重要的幕後功臣

1. pH值緩衝範圍: 6.8-8.2 · 涵蓋細胞的最佳培養環境: pH7.2-pH7.4
2. pH受溫度影響小 · 可貫串整個細胞及後續萃取實驗
3. 低溫時也能保持酶的結構和功能活性

限時九折優惠中(~110/12/31)

隨貨附COA, 純度>99.9%

100g / ~~980元~~ 882元

500g / ~~3500元~~ 3,150元



CAS NO.	品項	包裝	價格
4432-31-9	MES	100g/ 500g	2,000/ 7,560
145224-94-8	MES-H	100g/ 500g	2,000/ 7,560
1132-61-2	MOPS	100g/ 500g	1,050/ 3,675
5625-37-6	PIPES	100g/ 500g	2,100/ 7,350
29915-38-6	TAPS	100g	1,680
7365-44-8	TES	100g/ 500g	1,890/ 7,560
77-86-1	Tris Base	1Kg/ 5Kg	1,980/ 8,950
1185-53-141	Tris HCl	500g/ 500gx 2	3,360/ 5,250
5704-04-1	Tricine	100g/ 500g	1,050/ 4,200

預配滅菌緩衝劑溶液

2550	0.2M Bis-Tris Buffer, pH6.5, Sterile	100mL/250mL	1,050/1,880
------	--	-------------	-------------



Cyrusbio 塞魯士生技

0800-007-101/02-22988180

Line ID: @cyrusbio



最方便
的購買
管道

分生實驗預配緩衝溶液 (綠色字體產品為Western Blot專用)

型號	品項	包裝	價格
2003-50	<u>TAE Buffer 50X</u>	1L/ 4L	1,100/3,350
2002-5	<u>TBE Buffer 5X</u>	1L/ 4L	600/1,750
2007-10	<u>PBS Buffer 10X, with 0.1% preservative</u>	1L/ 4L	500/1,600
2007-10-S	<u>PBS Buffer 10X, Sterile</u>	1L	1,260
2538-10	<u>DPBS Buffer 10X, w/ calcium and magnesium, Sterile</u>	500mL	1,450
2539-10	<u>DPBS Buffer 10X, w/o calcium and magnesium, Sterile</u>	500mL	1,250
20074	<u>TBS Buffer 10X</u>	1L/ 4L	600/1,750
2532	<u>TBST Buffer 10X, Sterile (1X with 0.05% tween 20)</u>	500mL/ 500mLx2	1,100/1,680
2533	<u>TBST Buffer 1X, Sterile (with 0.05% tween 20)</u>	500mL/ 500mLx2	980/1,450
2512	<u>Upper Gel Buffer 4X</u>	1L	1,350
2522	<u>Lower Gel Buffer 4X</u>	1L	1,350
200A001	<u>TG Buffer (Transfer Buffer) 10X</u>	1L/ 4L	600/1,750
20080	<u>TGS Buffer (Running Buffer) 10X</u>	1L/ 4L	600/1,750
2501	<u>TE Buffer 10X</u>	1L/ 4L	600/1,750
2503	<u>0.5M Tris-HCl Buffer, pH6.8</u>	1L	800
2504	<u>1.5M Tris-HCl Buffer, pH8.8</u>	1L	800
2506	<u>1M Tris-HCl Buffer, pH6.8</u>	1L	1,050
2528	<u>3M Sodium Acetate, pH5.2, Sterile</u>	250mL/ 500mL	1,200/1,800
2502	<u>0.5M EDTA, pH8.0, Sterile</u>	1L	1,200
2511	<u>10% SDS Solution</u>	250mL/ 500mL	700/1,200
2514	<u>5M Sodium Hydroxide</u>	100mL	980
2530	<u>10% NP-40 Substitute Solution, Sterile</u>	500mL/ 500mLx2	1,200/1,800
2512	<u>1M Glycine, pH2.8, Sterile</u>	500mL/ 500mLx2	1,200/1,800
2531	<u>2M Magnesium chloride solution, Sterile</u>	500mL/ 500mLx2	800/1,350
2536	<u>20% PEG 6000, in 2.5 M NaCl, Sterile</u>	500mL/ 500mLx2	980/1,580



Cyrusbio 塞魯士生技

0800-007-101/02-22988180

Line ID: @cyrusbio



最方便的
購買
管道